



Safe Place Bedding LLC

Safe Place Travel Bed
Basic Set-up

Folding Your
Safe Place Travel Bed

Air Chamber Replacement
for Safe Place Travel Bed



INSTRUCTIONS



Safe Place
Bedding LLC

Setting Up Your Safe Place Bed

1

Each bundle includes a Safe Place Travel Bed, instructions, 4 straps, a duffel bag, an electric air pump and a custom machine-washable incontinence pad.




2

Place travel bed on full (74" W x 54" L) or larger stationary bed. Unfold bed long ways on top of stationary bed. Ensure the netting on the door closest to you is fully visible.



3

Push zipper head up out of zipper pocket on both sides of netting. Open zipper until zipper is 1" above valve. (Do not open zipper more than 1 inch above valve. Doing so will make zipper difficult to close after inflation.)




4

Ensure the black valve is completely inserted into the clear opening with no gap between the top of the clear plastic and the black valve. Then, repeat step 3 and 4 for the inflation valves located at the base, head, and bilateral long sides of the bed.



5

Place the tip of air pump into black valve at inflation point. You may need to hold tip of air pump in black valve during inflation until wrinkles are out of plastic. Repeat this process on all four chambers. All four chambers should be inflated and bed should be erected.




6

To achieve the desired look, place the inflated bed in the center of the stationary bed with equal distance on both sides from the edges of the bed. There should be at least 2-3 inches between the edge of the travel bed and side of mattress. Using bed on smaller than full mattress can lead to serious injury or death. Refer to the picture for guidance.




7

Insert the clear plug into the black one-way valves, then press down on the valve assembly into the inflatable chamber. Repeat this process for all four valves. The valve should be flush with the sides of the inflatable chamber.



Zip chamber closed and push zipper head back into pocket.



INSTRUCTIONS



Safe Place
Bedding LLC

Setting Up Your Safe Place Bed

8

To attach straps to your bed, connect first strap to the long side of the travel bed located directly below netting door on long side, loop strap under bed frame to opposite side of bed and snap into place. Repeat for second strap to long side of the bed. Next attach third strap to small triangle side of bed, loop strap longways under bed frame to opposite small triangle side and attach strap. Repeat for fourth strap. Follow step 10 to tighten straps attaching bed to frame. Not tightening straps can result in serious injury. During every use, please inspect strapping and buckles for damage.



9

To lock the bed zipper, connect the latch to the zipper. When inflating and deflating the bed, please ensure the doors are closed and locked.



Daily Checks While Using Travel Bed:

- Before placing user in bed, do a visual inspection of any imperfections in netting (holes, etc.). Stop using bed and contact customer service if netting has any imperfections.
- Manually inspect by squeezing air chambers of bed (4 total) to ensure they are still inflated. If one chamber is deflated, take out clear plug and black one-way valve from clear plastic opening/inflation point and follow steps 4-7 of instruction manual. Do not use bed for 24 hours to ensure chamber will hold air. If the chamber deflates again, contact customer service.
- Inspect strapping by pulling on each strap with your hand to ensure each strap is tight and bed remains in center of full or larger bed.
- Inspect zippers as you close the door. If you see a problem, do not use the bed and contact customer service. Check zipper lock before using the bed each night to ensure zipper lock has not loosened. To ensure the zipper lock is secure, no gap should be seen between fabric and zipper.
- Clean bed with warm water and a mild detergent on cloth if needed. We highly recommend using an incontinence pad to help keep bed clean. The included incontinence pad is machine washable. Contact customer service to purchase additional incontinence pads.

10

Tighten Bed Straps:

Work from side to side of the bed removing slack from each strap until the bed is firmly attached to the stationary bed.

***Note that there are two buckles on each strap.**

To tighten zipper lock and bed straps, pull strapping up in the middle of the buckle taking slack away above the buckle. To lock strapping in place, pull the tail of the strapping.

To ensure the zipper lock is tightly holding the zipper in place, repeat process until there is no gap between zipper and cloth.



Changing an Inflatable Chamber:

To change triangle-end chamber:

1. Open zippers at both sides of netting triangle.



2. Pull end of triangle chamber out of chamber cavity.



3. Pull inflatable chamber out on opposite side of triangle.



4. Reverse process to reinstall chamber.

To Change Long Chamber:

1. Unzip both zippers on bilateral sides of netting. Pull end of long chamber out of chamber cavity.



2. Push inflatable chamber into top canvas cavity.



3. Move to the opposite side and pull inflatable chamber out of opposite side.



4. Reverse process to reinstall inflatable chamber.

When reinstalling chamber make sure valve is pointed outward during reinstalling chamber.

WARNINGS

Stay Safe While Using Your Bed



**Safe Place
Bedding LLC**

Risk of suffocation if used before 36 months old or low muscle tone. Please consult your Doctor before using production.



Warning: Not suitable for children under 36 months. Small parts.



No smoking or heat source around products, flammability has not been tested.

- Always consult your medical Doctor before using this product.
- Keep away from water, not a flotation device.
- If one side of bed deflates, stop using bed.
- Firmly attach to stationary bed when in use, not tested independently.
- Sleeping mat is not a tested accessory, speak with your medical doctor and physical therapist before use.
- Sale Place Bed is twin (39" W x 75" L) in size, must be used on a full-size bed (74" W x 54" L) or larger bed.
- Zipper lock must be tightened before use and then checked nightly to prevent entrapment. If the zipper lock strap fails to remain tight, stop using bed, and contact Safe Place Bedding.
- Monitor user closely when he/she has access to outside of bed, risk may include but not limited to strangulation from straps, choking on zipper and entrapment.

Keep sharp objects away from product, punchers will not be covered in warranty and may cause injury to child.

