EasyStand Glider Moderate Support Package Feature Description

FEATURE	IMAGE	DESCRIPTION
Adjustable Foot Plates		Foot Plates are independently adjustable in height.
Secure Foot Straps		Velcro® and D-ring adjustment hold feet securely in place at ankle and toes. Must choose size in drop box.

FEATURE	IMAGE	DESCRIPTION
Flip-Up Knee Pads	Flip-Up Knee Pads	Knee Pads Flip-Up and out of the way for transfers. They are adjustable in height.
Full Range Seat	Fill Range Seat with hinged Field-away sections for each	The Full Range Seat has hinged, break-away sections for each leg that allow full hip and knee extension with ankle dorsiflexion. The seat is wide to enhance transfers.

FEATURE	IMAGE	DESCRIPTION
Easy Adjust Seat Depth	Eas/Stand	Roller bearing seat depth allows for smoother seat depth adjustment.
Velcro [®] Positioning Belt		2" (5cm) wide positioning belt has a D-ring Velcro® closure, provides hip stability for user. Fits hip circumference of 26"-55" (66- 140cm) OR 38"-67" (97-170cm). Must choose size in drop box
Hip Supports		Supports the hips in symmetrical alignment from sitting to standing. Pad size is 4"x6" (10x16cm). Must choose width range of hip supports in drop box.
Removable Flat Back		Removable back releases with two knobs and provides rear access for patient lift transfers. Must choose height of back in drop box.

FEATURE	IMAGE	DESCRIPTION
Chest Strap		Velcro® and D-ring adjustment. Chest strap pad is 5"Wx19"L (13x48cm). Must choose size of strap in drop box.
Tray and Chest Pad		The Clear Tray and Chest Pad are adjustable in depth and height.
Adjustable/ Removable Actuator Handle		Removable actuator handle adjusts the distance of reach to actuator handle to meet user's needs.

FEATURE	IMAGE	DESCRIPTION
Handle Extensions- For Medium Gliders Only		U-shaped extensions bring handles 6" (15cm) closer to the user for positioning and range. Handle height is adjustable. For Medium Glider only.
Handles for Large Gliders Only		Handle Height is easily adjustable. For Large Glider only.
Adjustable Resistance Cylinders	<image/>	The Adjustable Resistance Cylinders are accessible to the user and have 12 different settings. It allows the user to start with minimal resistance and increase as they become stronger or they can be unhooked completely.

FEATURE	IMAGE	DESCRIPTION
Positioning Bar		The Positioning Grab Bar provides a place for a sturdy grip during transfers.
5" Front Wheels and 5" Rear Locking Casters	5" Rear Locking Casters	Front wheels are fixed and roll forward and back. Back wheels are locking swivel casters.