

Grillo front drive sizes and weight

- A: height ground to thoracic support
- B: height ground to pelvic support
- C: length encumbrance
- D: overall width
- E: width of inner supports
- F: circumference of supports
- G: user height
- Total weight
- Maximum load

mini PTM

from 60 to 76 cm
 from 42 to 61 cm
 74,5 cm
 62 cm
 from 19 to 25 cm
 from 53 to 79 cm
 from 75 to 95 cm
 15,4 kg
 35 kg

mini

from 51 to 75 cm
 from 42 to 61 cm
 74,5 cm
 62 cm
 from 19 to 25 cm
 from 53 to 79 cm
 from 75 to 95 cm
 15 kg
 35 kg

small




from 66 to 100 cm
 from 49 to 84 cm
 75 cm
 69 cm
 from 19 to 25 cm
 from 60 to 86 cm
 from 90 to 135 cm
 19 kg
 45 kg

medium

from 73 to 124 cm
 from 54 to 100 cm
 91 cm
 76 cm
 from 22 to 36 cm
 from 74 to 120 cm
 from 110 to 165 cm
 22,5 kg
 90 kg

large

da 84 a 147 cm
 da 63 a 120 cm
 111 cm
 80 cm
 da 24 a 40 cm
 da 81 a 140 cm
 da 150 a 190 cm
 29,5 kg
 110 kg



Make sure that measures A and B indicated in the figure do not exceed the suggested range indicated for the supports.

Grillo rear drive sizes and weight

B: height ground to pelvic support
 C: length encumbrance
 D: overall width
 E: width of inner supports
 F: circumference of supports
 G: user height
 Total weight
 Maximum load

mini

from 42 to 61 cm
 74,5 cm
 62 cm
 from 19 to 25 cm
 from 53 to 79 cm
 from 75 to 95 cm
 15 kg
 35 kg

small



from 49 to 84 cm
 75 cm
 69 cm
 from 19 to 25 cm
 from 60 to 86 cm
 from 90 to 135 cm
 18,5 kg
 45 kg

medium

from 54 to 100 cm
 91 cm
 76 cm
 from 22 to 36 cm
 from 74 to 120 cm
 from 110 to 165 cm
 21,7 kg
 90 kg

large

from 63 to 120 cm
 111 cm
 80 cm
 from 24 to 40 cm
 from 81 to 140 cm
 from 150 to 190 cm
 28,5 kg
 110 kg



Make sure that measure B indicated in the figure do not exceed the suggested range indicated for the pelvic support.